

PULSAR SCIENTIFIC

Your physician has recommended that you utilize a Compression Accelerated Recovery System following your upcoming surgery. The system is designed to provide exceptional post-operative recovery and pain relief which has been shown to translate into less narcotic use to control pain.

Why Do COMPRESSION & COLD Work So Well Together?

Studies show that cold therapy decreases pain, muscle spasms, and swelling in soft tissue injuries and decreases the likelihood of tissue damage. Compression pushes excess fluid away from the injury, which helps to eliminate existing swelling and minimize new swelling. When combined, these two therapies produce a new level of healing.

Recovery+ delivers a unique combination of compression and cooling, as well as comfortable, form-fitting anatomic wraps, and adjustable pressure and temperature settings, all in a state-of-the-art system. Simply put, there is nothing else like it.



8 Stony Brook Lane
Ludlow, MA 01056

Phone: 413.589.0851
Fax: 413.385.0003



INTRODUCING *Recovery+*™

PORTABLE DVT PROPHYLAXIS

6 COMPRESSION SETTINGS

CONSTANT COOLING

CONTRAST THERAPY



TREATMENT MODALITY AND INDICATIONS

TREATMENT MODALITY

Compression Therapy

6 Compression settings: 15 - 75 mmHg.

Thermal Therapy

43° to 50° temperature.

Portable DVT Prophylaxis

Decrease the risk of deep venous thrombosis (DVT).

INDICATIONS

Edema

Reduction of edema associated with soft tissue injuries such as burns, postoperative edema, and ligament sprains.

DVT Prophylaxis


DVTherapy offers greater range of compression, allowing physicians more flexibility in prescribing treatments.

Post Operative Orthopedic Surgery

Treatment of pain and swelling of acute periarticular processes, treatment of pain and swelling following mobilization of shoulder stiffness under anesthesia, treatment of pain and swelling postoperatively for bones, joints and soft tissue, treatment of pain and swelling caused by musculoskeletal contusions and athletic injury.

Pain

Cold therapy for pain management.



BREADTH OF TREATMENT



Therapies:

- Constant Cooling
- 6 Compression Settings
- Contrast Therapy
- Portable DVT Prophylaxis

WRAPS TO COVER ANY SURGERY

- | | |
|------------|------------|
| • Shoulder | • Knee |
| • Ankle | • Back |
| • Hip | • Elbow |
| • DVT Calf | • DVT Foot |

FREQUENTLY ASKED QUESTIONS

How long do I keep the Recovery+?

The doctor typically orders the Recovery+ for 21 days. If the unit is needed for a longer period of time, contact your doctor's office and ask that your prescription be extended.

Who do I call if I have questions?

Please call your local distributor or Pulsar Scientific, **413.385.0851**.

How do I diagnose the Recovery+?

The reliable design of Recovery+ provides easy-to-access, serviceable wear parts. The device is able to self diagnose when mechanical faults or operational mistakes are made, minimizing service calls and allowing DIY maintenance.

How often do I use Recovery+?

Please refer to the instructions provided by your physician.


Recovery+™

Optional Items:

- Bed Hanging Basket
- Dual Connector Hose

HEAL FASTER AND REDUCE PAIN WITH SEGMENTAL COMPRESSION.

1. Continuous water circulation and segmental compression promote healing.
2. Segmental compression (see arrows) forces blood flow away from surgical site to reduce swelling.
3. 100% latex-free wraps reduce pain, swelling, edema and hemarthrosis, and promote faster healing.



WMI ENTERPRISES, LLC

7750 Zionsville Road, Suite 850 B
Indianapolis, IN 46268


Toll Free: 800.592.2274

Phone: 317.704.3300

Fax: 317.704.3303

© 2013 WMI Enterprises, LLC. All wraps have patents granted or pending

INTRODUCING THE 2000 SERIES SEGMENTAL THERMAL COMPRESSION WRAPS



EXPERIENCE PRECISE COMPRESSION
& TEMPERATURE CONTROL

The Gentle Compression and Adjustable Cold
MEAN GREATER COMFORT AND BETTER OUTCOMES.

2000 SERIES A FULL LINE OF ERGONOMIC WRAPS

